

## Tuttinauer Rotary Sealing Machine

### TSM 850 Rotary Sealing Machine

The sterilization of medical instruments begins with the packaging. That is to say the packaging, and therefore also the validation of the packaging process, is crucial for guaranteeing sterility up to the point of use and for meeting quality standards. An appropriate packaging system facilitates the safe packaging, sterilization, storage and aseptic removal of the sterile goods. Sealing a sterile package such as this is only possible with a sealing device whose functions meet the requirements of the German Industrial Standard DIN 58953-7 and which can be validated in terms of ISO 9001:2000.

- Continuous monitoring of the relevant parameters: the unit is switched off automatically when there is a deviation from the set value
- Seal-check function: direct connection to the ht 150 SCD seal checking device

### Microprocessor-controlled printing unit

Printing takes place during the sealing process over the sealed seam:

- Sterilization and use by date with automatic updating
- Batch and personnel number
- Symbols in accordance with EN 980
- CE mark with subsequent text (for medical products in accordance with the law on medical products)
- Text memory
- Variable font size
- Printer can be switched off



### Functions

- Menu-based operation via digital LCD display
- Standard texts in German, English, French, Italian and Spanish
- Clock and calendar function
- Data storage and date extender even when device is switched off
- Automatic maintenance check
- "Seal-Check" function for the continuous process control
- Follow-up calibration service with certificate
- Standby function: the device goes into energy-saving standby mode after a predefined time (10 - 120 min)
- Batch counter can be reset

## Technical data TSM 850

Printing unit	single line
Print start from edge	infinitely adjustable in mm
Sealing speed	10 m/min
Sealing edge,	infinitely adjustable 5 – 30 mm
Sealed seam width	12 mm
Safe dist. betw. seam & sealed object	(DIN 58953-7) (A) > 30 mm
Temperature control	microprocessor
Sealing temperature	80 – 220°C
Temperature tolerance	< ± 1 %
Power	500 VA
Mains connection	230 V/50 Hz (115 V/60 Hz)
Size w x d x h	620 x 260 x 250 mm
Weight	23 kg
Housing	chrome steel

## Accessories

- **Roller table TSM 500 RT**, dimensions (wxhxd) 870 x 310 x 40mm
- **TSM 611 S Cutting device with roll holder in chrome steel**

Dimensions hm 611 S (wxhxd) 458 x 125 x 344 mm  
Roll diameter max. 200 mm  
Cutting length max. 400 mm

